

# Quick Guide Advanced Modes **C-FLOW / MPP & MPV\***



\* C-FLOW: Continuous Flow Mode  
MPP: Mouthpiece Ventilation  
by Pressure  
MPV: Mouthpiece Ventilation  
by Volume

This guide does not  
replace the User Manual

Intended for all the people  
authorized to adjust settings

# C-FLOW (Continuous Flow) Mode

## PRINCIPLE:

It is a high-flow humidified therapy that provides respiratory support to spontaneously breathing patients by delivering heated and humidified air at a flow rate up to 60 LPM.

Supplementary oxygen is possible in addition to the therapy.


It is possible to set the Pmax, if reached will not exceed the set limit.

## Setup the C-FLOW mode



Select **CLINICAL**, then on the tab **CONFIG**, select:

- 1 **Circuit type:** Leak
- 2 **Ventilation:** with **Flow** type and **C-Flow** mode.



### 3 Set the parameters on the tab **SETTINGS**:

Settings	Adult	Pediatric
Flow	10 - 60 L/min	2 - 25 L/min




## Supplemental oxygen - % FiO<sub>2</sub>

FiO <sub>2</sub> % theoretical value	C-Flow setting (L/min)							
		10	15	20	30	40	50	60
O <sub>2</sub> Flow (L/min)	4	53	42	37	32	29	27	26
	8	84	63	53	42	37	34	32
	12	-	84	68	53	45	40	37
	16	-	-	84	63	53	46	42
	20	-	-	-	74	61	53	47
	25	-	-	-	87	70	61	54
	30	-	-	-	-	80	68	61
	40	-	-	-	-	-	84	74
	50	-	-	-	-	-	-	87

The oxygen flow must not exceed the set flow rate.

## Circuit configuration (only with a single-limb circuit)




How to set up from the humidifier:

- 1 Connect the short circuit to the inspiratory port of the ventilator.
- 2 Connect the heat circuit, temperature probes and electric adapter.
- 3 Connect the nasal cannula.

\* Use the leak adapter (cap) included into the starter kit.



- 4 Select the most appropriate setting on the humidifier.
- 5 Fit the nasal cannula to the patient



## Compatible Accessories and Consumables

		Hospital		Homecare	
		Fisher & Paykel HEALTHCARE	WILAMED	Fisher & Paykel HEALTHCARE	WILAMED
Adult	Humidifier model	MR850	AIRcon	HC550	AIRcon
	Circuit kit (including chamber)	RT202	BTS3234A	RT202	BTS3234A
	Nasal cannula	Optiflow™	WILAhiflow	Optiflow™	WILAhiflow
Pediatric	Humidifier	MR850		HC550	
	Circuit kit (including chamber)	RT331		RT331	
	Nasal cannula	Optiflow™		Optiflow™	

**WARNING:** The circuit diameter 22 mm (between the ventilator and the humidifier) is not included in the pediatric kit.


## MPP & MPV

Mouthpiece Ventilation is a form of non-invasive ventilation via a mouthpiece, dispensing the need for face or oral masks. It is indicated for patients who require ventilatory support and have good cooperation levels.

### Mode MPP : Mouthpiece Ventilation in Pressure

Frequency = OFF  
 Min. pressure = 5 mbar  
 Pres. Control = 25 mbar

- 1 = Patient-triggered cycle
- 2 = Disconnection alarm



Settings	Adult	Pediatric	Limits
Pres. control (mb)	5-49	5-49	-
Bias Flow (l/min)	OFF / 1-20 l/min	OFF / 1-20 l/min	-
Slope	1-5 (100-500 ms)	1-5 (50-250 ms)	-
Frequency (cpm)	OFF / 5-60	OFF / 5-80	Frequency $\leq 30$ / Insp. time (I/E $\leq 1/1$ )
Inspiratory time (s)	0,3-2,5	0,3-2,5	Frequency $\leq 30$ / Insp. time (I/E $\leq 1/1$ )
I. Trig	AUTO / 1-5	AUTO / 1-5	-

### Articulated instruments reference:

Distributor	Reference
Lowenstein	WM27647
Philips	1102862


### Mouthpiece reference:

Distributor	Reference
Lowenstein	WM27646

## MPV Mode : Mouthpiece Ventilation in Volume

Volume = 600 ml  
 Rate = OFF  
 Min. pressure = 5 mbar

- 1 = Patient-triggered cycle
- 2 = Disconnection alarm



Settings	Adult	Pediatric	Limits
VT (ml)	100-2500	100-600	-
Blas Flow (l/min)	OFF / 1-20 l/min	OFF / 1-20 l/min	-
Flow ramp	1.Square, 2.Decelerating		-
Frequency (cpm)	OFF / 5-60	OFF / 5-80	Frequency $\leq 30$ / Insp. time (I/E $\leq 1/1$ )
Inspiratory time (s)	0,3-2,5	0,3-2,5	Frequency $\leq 30$ / Insp. time (I/E $\leq 1/1$ )
I. Trig	AUTO / 1-5	AUTO / 1-5	-

### Manufactured by

EOVE  
 Helioparc  
 4, Boulevard Lucien Favre  
 64000 Pau  
 France  
[www.eove.fr](http://www.eove.fr)



EO-150 provides continuous or intermittent ventilation support for adults and children weighing at least 3.5 kg needing mechanical ventilation, at home, in a facility or hospital, and in portable environments, for invasive or non-invasive ventilation.

Medical device class IIb - CE 0459 - Manufactured by EOVE. Please read the user manual.