

# Respireo

Primo F



## Simple & convenient facial mask



We are all different. The Respireo Primo series is a range of masks for ventilation that has been specifically designed to fit each need of the patients and adapt to each morphology.

Simple & convenient, Respireo Primo F is a full face mask with adjustment dial and high ease of use. It comes Vented and Non Vented and has been specially conceived to ensure your patients' maximum comfort. It is also recommended for patients with high therapy pressures, mouth breathing or with difficulty to handle things, thanks to the intuitive fastening clip and a simple assembly.




# Respireo

## Primo F



### Dual wall cushion

- Ensuring stability and airtightness with a great comfort thanks to the thin external lip
- Easy assembly with only 2 components




### Well distributed and adjustable forehead support



- Adjusting dial to find the perfect fitting with the forehead inclination
- Large silicone pad to distribute the pressure points

### "V" shaped clips

- Flared clip system for an intuitive fitting and a quick release system if need be



### Diffuse and silent vent



Respireo Primo F is also available in **non vented** version for patients with a valve ventilation.



Codes	Designation	Size
KM217900	Respireo Primo F (manufacturer reference 472530*)	S
KM218000	Respireo Primo F (manufacturer reference 472531*)	M
KM218100	Respireo Primo F (manufacturer reference 472532*)	L
KM218200	Respireo Primo F Non Vented (manufacturer reference 472533*)	S
KM218300	Respireo Primo F Non Vented (manufacturer reference 472534*)	M
KM218400	Respireo Primo F Non Vented (manufacturer reference 472535*)	L
KM218500	Cushion Respireo Primo F (x10) (manufacturer reference 472615)	S
KM218600	Cushion Respireo Primo F (x10) (manufacturer reference 472616)	M
KM218700	Cushion Respireo Primo F (x10) (manufacturer reference 472617)	L
KM224000	Headgear Respireo Standard (manufacturer reference 472456**)	Stand
KM224100	Headgear Respireo (manufacturer reference 472457**)	L
KM220400	Kit 20 sizing gauge Respireo Primo F (manufacturer reference 472432)	



Please read user manual.

(\*) Medical device class IIa CE 0051. (\*\*) Medical device classe I CE  
 manufactured by Air Liquide Medical Systems S.r.l.  
 Distributed by Air Liquide Medical Systems S.A.